

USER'S MANUAL

FOR OPERATION PRINTER MODEL CSP-160 II

SHINKO DENSHI CO., LTD.
Tokyo

3325M81

CONTENTS

1. INTRODUCTION.....	1
2. GENERAL SPECIFICATIONS	2
3. NAMES OF PARTS	3
4. The NAME AND FUNCTION OF ALL PARTS.....	3
5. INSTALLATION.....	5
6. ERROR MESSAGES	6
7. DATA OPERATIONS	7
8. OPERATIONS	8
9. CONNECTION WITH SCALE.....	10
10. CLOCK ADJUSTMENT.....	13

USRES'MANUAL
MODEL CSP-160 II

1. INTRODUCTION

The Operation Micro-Printer CSP-160II, storing data from your Shinko electronic scale, operates Total, Average, Deviation, Maximum, Minimum, Ratio of Data of out of Limits, etc., besides ordinary data recording.

The results are printed by cartridge ink ribbon.

1-1. Features

- * Statistical Operation and Printing

Number of Data	(N)	Average	(AVE)
Maximum	(MAX)	Standard Deviation	(σ_n , σ_{n-1})
Minimum	(MIN)	Failure Ratio	(P)
Range	(R)	Total	(TOTAL)
- * Histogram
- * Judgements: HIGH, LOW or GOOD of data
- * AUTO/MANUAL selectable

CAUTIONS ON USE

1.	Always use the attached AC adaptor only, by connecting to an isolated line outlet from other devices to avoid disturbing noises.
2.	Turn off the power switch of CSP-160II first, before cutting line power off. Never unplug the AC adaptor from the outlet before turning the printer switch off. Turn off power switches of both printer and scale before connecting/disconnecting cable and the AC adaptor.
3.	Do not print without setting roll paper and/or ink ribbon.
4.	Position and use the printer in a suitable place. A quiet area without vibration, dusts, humidity, heat or cooling sources, is preferable.
5.	After operation, turn off the power switch of the printer, then unplug the AC adaptor from The line outlet without fail.


2. GENERAL SPECIFICATIONS

Character Construction	: 5×7 Dot Matrices
Printing Digits	: 16 digits maximum
Printing Speed	: Approx. 1 line/sec.
Character Sizes	: 1.8mm(W)×2.5mm(H)
Operation Temperature	: 0°C ~ 40°C
Dimensions	: 155mm (W) ×145mm (D) ×48mm (H)
Power Source	: AC adaptor
Weight	: Approx. 400g without AC adaptor and roll paper
Ink Ribbon / Color	: Cartridge ERC-05 Violet
Recording Paper	: 45mm(W)×30mm diameter ordinary paper

The following items are included in the package. If there are any missing or damaged items, please contact our sales office or the dealer of purchase.


1. CSP-160II (Ink ribbon included) 1
2. Recording paper 2
3. Connecting cable Din5P 1
4. Connecting cable Dsub9P 1
5. Foot switch jack 1
6. Instruction manual 1

3. NAMES OF PARTS




4. The name and function of all parts

4.1 Front part



①	POWER lamp	Lights on to indicate power is on.	
		Also, indicates by colors:	
		Green:	In normal condition of printer.
		Red:	Error in printer mechanism, such as locked by paper jam.
	Orange:	Error in internal clock.	
②	DATA	Key to input data manually.	
③	FEED	Key to feed recording paper in inching by one line.	
④	STAT	Key to process statistical data, printing results and histogram, by the procedure described in page 8 to page 9.	
⑤	LIMIT	Key to start/finish setting of limits.	
⑥	CANCEL	Key to cancel the latest input data by printing "CANCEL".	
⑦	RESET	Key to clear all data including LMIT values and operation results.	
④	STAT	Key to process statistical data, printing results and histogram, by the procedure described in page 8 to page 9.	

4.2 Rear side, Right side



① POWER SWITCH	Switch for on/off of the power to the CSP-160II
② S.IN jack	Jack for data from scale by connecting the attached cable. Corresponding plug: T-0718 by TEISHIN
③ F.SW jack	Jack for foot Switch to operate the same performance as DATA key. Corresponding plug: T-0706 by TEISHIN
④ Dip Switches	Switches for setting function of the printer.
⑤ DC IN JACK	Jack for connection of attached AC adaptor.

SETTING OF DIP SWITCHES

SW №	ON	OFF	Setting at the time of delivery
1	Baud Rate Setting, See the table below.		OFF
2			OFF
3	Command by scale	Command by printer	ON
4	Statistical Calculation	-	OFF
5	English character	Japanese character	ON
6	Totaling	Statistical Calculation	ON
7	Automatic printing	Manual Printing	OFF
8	5 scales print reset	50 scales print reset	OFF

SW 1	SW 2	BAUD RATE
OFF	OFF	1200
ON	OFF	2400
OFF	ON	4800
ON	OFF	9600

- CAUTION** : 1. Turn power switch OFF before setting.
 2. SW No.9 and 10 are not in use.
 3. SW No.4 ON to use the statistical calculation function when the scale is displaying the auxiliary scale interval


5. INSTALLATION

5.1 Connections


1. Verify that power switches of the CSP-160II Printer and the scale are both OFF.
2. Connection with the scale. Connect one of the connectors of the attached cable with S.IN jack of the CSP-160II and the other connector of the cable to output jack of the scale.
3. Connection of the AC adaptor.
Connect the jack of attached AC adaptor with DC.IN jack of the CSP-160II.
4. Connection of foot switch option
Connect the foot switch F.SW jack of the CSP-160II.

5.2 Ink Ribbon Setting

1. Remove the printer cover by lifting it up.
2. To remove old cartridge, press the end of the cartridge at "PUSH" mark.





3. Setting new ribbon to platen, put a new cartridge in, fitting knob axis. Turn the knob to arrow direction to get rid of the slack in the ribbon.



NOTE : The ink ribbon is disposable, don't put ink to the ribbon.

5.3 Roll paper Setting

1. Under power ON, remove the printer cover by lifting it up.
2. Cut the paper tip straight.
3. Inserting the tip into the paper slot, press the "FEED" key to feed the paper out of the outlet.
4. Position the roll paper in the holder, then cover the printer after feeding the paper out of the opening of the Cover.



6. ERROR MESSAGES

- | | |
|-------------|--|
| LIMIT ERROR | : Upper and lower limit are set at the same value, or set only one limit. |
| I/O ERROR | : Incorrect connection of the printer with scale, or data is out of specification. |
| STAT ERROR | : STAT key is pressed with only one data inputted. |
| DATA ERROR | : Data from the scale is out of range. |
| ERROR | : Printer malfunction by disturbance noise, etc. |

7. DATA OPERATIONS

7.1

ITEMS TO BE PRINTED	SETTING CONDITION OF STATISTICAL OPERATION		
	DIP SWITCH NO.6 SET "ON"	DIP SWITCH NO.6 SET "SET OFF"	
		No Limit are set	Limits are set
TOTAL	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
N: Number of data	-	<input type="radio"/>	<input type="radio"/>
MAX: Maximum	-	<input type="radio"/>	<input type="radio"/>
MIN: Minimum	-	<input type="radio"/>	<input type="radio"/>
R: Range	-	<input type="radio"/>	<input type="radio"/>
AVE: Average	-	<input type="radio"/>	<input type="radio"/>
Standard deviation	-	<input type="radio"/>	<input type="radio"/>
P: Failure Ratio	-	-	<input type="radio"/>
Histogram	-	-	<input type="radio"/> *1
Judgement result	<input type="radio"/> *2	-	<input type="radio"/>
Date	-	<input type="radio"/>	<input type="radio"/>

"○" mark shows items to be printed.

*1: Effective when limits are set by the CSP-160II Printer, not by scale.

*2: Effective when limits are set.

7.2 OPERATION FORMULAE

Total (TOTAL) = $\sum x_i$

Range (R) = (MAX) - (MIN)

Average (AVE) = $\sum x_i / N$

Standard Deviation (σ_n) = $\sqrt{\{N \cdot \sum x_i^2 - (\sum x_i)^2\} / N^2}$

Standard Deviation (σ_{n-1}) = $\sqrt{\{N \cdot \sum x_i^2 - (\sum x_i)^2\} / N \cdot (N - 1)}$

Failure ratio (P) = (HIGH) + (LOW) / N

7.3 JUDGMENT FORMULAE

Upper limit Value < Data + NG (Prints "H")

Upper limit \geq Data \geq Lower Limit GOOD (space)

Upper limit Value > Data - NG (Prints "L")

8. OPERATIONS

PRINT SAMPLE

```

*RESET
*LIMIT DATA
LMT1      599.7 g
LMT2      2199.4 g
*SET LIMIT DATA
HILM      2199.4 g
LOLM      599.7 g

001      1199.8 g
002      999.8 g
003      599.8 g
004      1999.5 g
005      1599.4 g
006      1199.5 g
H007      2399.4 g
008      1399.5 g
009      799.7 g
010      599.8 g
L011      399.8 g
012      1599.3 g
L013      1002.7 g
014      802.7 g
    
```

```

024      1199.8 g
025      999.8 g
026      1199.8 g
027      1199.9 g
028      1199.6 g
029      1200.7 g
030      1200.8 g
031      1200.9 g
032      1199.4 g
033      999.8 g
034*     2004.1 g
035      1199.7 g
036      1199.7 g
L037      199.9 g
*CANCEL
037      999.7 g
038      799.4 g
039      599.8 g
040      799.8 g
    
```

KEY OPERATION

RESET

.... To clear memories

LIMIT

.... To set limits
Press for judgements & histogram.

DATA

.... Input one limit value

DATA

.... Input another limit value

LIMIT

.... To finish setting of limits

DATA

.... Data printing, by either of Data key of Printer, or by foot switch, or by command from scale.

DATA

Prints * mark with an unstabilized data.


Mark	Judgement Result
L	Data is less than limit
H	Data is over than limit

CANCEL

.... To cancel the latest data

PRINT SAMPLE

KEY OPERATION



9. CONNECTION WITH SCALE


In case the printer CSP-160II is used in connection with any SHINKO electronic scale, set both the printer and the scale referring to the following. Also refer to operation manual for the scale.

9-1 For Automatic Printing

: Printing starts automatically by stabilization of data from the scale.

1. Setting of the Printer

Set dip switches on the rear of the CSP-160II as follows:



Dip Sw. No. 6

OFF : Statistical Operation

ON : Totaling only

Dip Sw. No.9 & 10 = OFF

2. Setting of the Scale Functions

Scale Series	Functions	Output Provision
AJ	7. IF. 1 71. O.C. 1	Standard equipment
HG, PF	8. IF. 1 81. O.C. 1	Standard equipment
CUXII	4. O.C. 1	Option
CGX, CGX-K	6. IF. 1 61. O.C. 1	Option
CG, CG-K, CG-KF	6. IF. 1 61. O.C. 1	Option
DJ-E	6. O.C. 1	Option
GZII	5. IF. 1 51. O.C. 1	Option
HTR, HT	7. IF. 2 71. O.C. 1	Standard equipment
LF	7. IF. 5 71. O.C. 1	Standard equipment
ALE, HTN	4 1 1 ACTIVATE: ON 4 1 2 FORMAT: CSP 7 4 1 3 CONDITION: 1	Standard equipment
AB	7. IF. 1 71. O.C. 1	Standard equipment
LN	6. IF. 1 61. O.C. 1	Standard equipment
CJ	6. IF. 1 61. O.C. 1	Standard equipment
SJ	7. IF. 1 71. O.C. 1	Option
HJ-K	6. IF. 1 61. O.C. 1	Standard equipment
FJK	6. IF. 1 61. O.C. 1	Standard equipment
DME	A. O.C. 1	Standard equipment
CT	7. IF. 1 71. O.C. 1	Standard equipment
FS, FZ	411. IC. 2 or 3 412. OO. 1 413. OC. 1 415. BL. 1 416. PA. 0	Standard equipment

9-2 For Manual Printing

A. Print command by DATA key of the CSP-160 II

Printing starts by every pressing of the DATA key, regardless of the status of data.

1. Setting of the Printer

Set dip switches on the rear of the CSP-160II as follows:

OFF		<u>Dip Sw. No. 6</u>
ON		OFF : Statistical Operation ON : Totaling only
		<u>Dip Sw. No.9 & 10 = OFF</u>


2. Setting of the Scale Functions: The same as above 9-1 2.

B. Print command by Scale Key

Printing starts by every pressing of a scale key, regardless of the status of data.

1. Setting of the Printer

Set dip switches on the rear of the CSP-160II as follows:



Dip Sw. No. 6

OFF : Statistical Operation

ON : Totaling only


Dip Sw. No.9 & 10 = OFF

2. Setting of the Scale Functions


Scale Series	Functions	Scale Key
AJ	7. IF. 1 71. O.C. 3	[Memory]
HG, PF	8. IF. 1 81. O.C. 3	[S]
CUXII	Impossible to set	
CGX, CGX-K		
CG, CG-K, CG-KF	6. IF. 1 61. O.C. 3	[S]
DJ-E	6. O.C. 3	[Print]
GZII	5. IF. 1 51. O.C. 3	[Print]
HTR, HT	7. IF. 2 71. O.C. 3	[Print]
LF	7. IF. 5 71. O.C. 3	[Print]
ALE, HTN	4 1 1 ACTIVATE: ON 4 1 2 FORMAT: CSP 7 4 1 3 CONDITION: 3	[Output]
AB	7. IF. 1 71. O.C. 3	[Print]
LN	6. IF. 1 61. O.C. 3	[Print]
CJ	6. IF. 1 61. O.C. 3	[Output] or [Print]
SJ	7. IF. 1 71. O.C. 3	[Memory]
HJ-K	6. IF. 1 61. O.C. 3	[Print]
FJK	6. IF. 1 61. O.C. 3	[Print]
DME	A. O.C. 3	[Print]
CT	7. IF. 1 71. O.C. 3	[Print]
FS, FZ	411. IC. 2 or 3 412. OO. 1 413. OC. 3 415. BL 1 416. PA 0	[Transfer]

10. CLOCK ADJUSTMENT

10-1 Setting Date and Time



10-2 Date and Time printing for checking



Press DATA key for more than 2 seconds

