

Specifications

EXTERNAL WEIGHT MODELS

| Model ※1 | LNA 623 | LNA 1202 | LNA 2202 | LNA 3202 | LNA 4202 | LNA 6202 | LNA 10002 | LNA 15001 | LNA 21001 | LNA 31001 |
|---------------------|---|----------|---------------|-----------|----------|----------|---------------|-----------|-----------|-----------|
| Capacity | 620g | 1200g | 2200g | 3200g | 4200g | 6200g | 10000g | 15000g | 21000g | 31000g |
| Read-out(d) | 0.001g | | | 0.01g | | | 0.05g | 0.1g | 0.1g | |
| Verification(e) ※2 | 0.01g | | | 0.1g | | | — | 1g | 1g | |
| Repeatability(s) | 0.001g | | | 0.01g | | | 0.05g | 0.1g | 0.1g | |
| Non-Linearity(typ.) | ±0.001g | | | ±0.01g | | | ±0.05g | ±0.1g | ±0.1g | |
| Pan size | 120×140mm | | | 200×200mm | | | | | 220×250mm | |
| Calibration | with external weight only | | | | | | | | | |
| Dimensions | 330×220×190mm (including windshield) | | 333×220×88mm | | | | 330×220×111mm | | | |
| Weights | approx. 3.5kg | | approx. 4.0kg | | | | approx. 8.5kg | | | |

INTERNAL WEIGHT MODELS

| Model ※1 | LNA 623R | LNA 1202R | LNA 2202R | LNA 3202R | LNA 4202R |
|---------------------|---|-----------|---------------|-----------|-----------|
| Capacity | 620g | 1200g | 2200g | 3200g | 4200g |
| Read-out(d) | 0.001g | | 0.01g | | |
| Verification(e) ※2 | 0.01g | | 0.1g | | |
| Repeatability(s) | 0.001g | | 0.01g | | |
| Non-Linearity(typ.) | ±0.001g | | ±0.01g | | |
| Pan size | 120×140mm | | 200×200mm | | |
| Calibration | with internal and external weight | | | | |
| Dimensions | 330×220×190mm (including windshield) | | 333×220×88mm | | |
| Weights | approx. 4.0kg | | approx. 6.0kg | | |

Options

| | |
|-------|------------------------------------|
| LNALM | Relay contact |
| LNAUH | Under weighing hook (1200g-15000g) |
| LNAR4 | RS422A output |
| LNADK | Density measurement kit |

Common Specification


Power source : AC120/230V, DC12V
 Output : RS232C (2 outputs)
 Measuring system : Tuning-fork frequency system
 Tare : Full weighing range
 Display : Fluorescent display

※1 followed by "CE" for Approval Type
 ※2 for Approval Type

What makes the tuning-fork sensor so precise?


The tuning-fork sensor measures force or mass by gauging changes in oscillation frequency when a load is applied to a long, narrow vibrator, and it digitally outputs the readings.

Unlike load cell or electromagnetic systems, the tuning-fork sensor does not rely on material distortion, electromagnetic force, heavy power consumption, or A/D converters, so its inherent margin of error is extremely small, and its high precision can be maintained for a long time.



Double-Ended tuning fork (DETF) vibrator

Dimensions



The contents of this catalogue are subject to change due to modifications and/or other reasons.

SHINKO DENSHI CO., LTD.

SHINKO DENSHI CO., LTD.
 1-52-1 ITABASHI, ITABASHI-KU, TOKYO
 173-0004, JAPAN
 TEL : +81-3-5944-1643 FAX : +81-3-6905-5526
 URL : http:// www.vibra.co.jp/global

Distributed by:



2308

LNA Series



Complete weighing solution with wide range and much applications

Statistic function for quality control use

Comfortable operation in quick response and stable indication

For laboratory, light&heavy industry, jewelry shops, etc...



PRECISION TUNING-FORK BALANCE

LNA Series



Statistic function for QC & etc...

$$Q_n = \frac{s}{\sqrt{n}} \cdot t(1-a)$$

ViBRA LNA series has the function to automatically calculate the various statistical data from the measurement results. This data can be useful for the quality control in the assembling line and the statistical checks of prepackaged products. The data can be output to the printer and PC.

Fluorescent display, clearly visible

The large fluorescent display is clearly visible. It can make it easy to operate the balance even in the dark location.



Quick response and stable indication

The quick response and the stable indication are important for almost all the weighing operations. ViBRA LNA series promises you the quickness and stableness so that it can make the measurement works much more efficient and less time-consuming.



Accurate measurement by appropriate calibration

It is highly important to to keep the accuracy of the balance by calibration. The procedure of the calibration is sometimes bothering, but in ViBRA LNA series, you can adjust the balance with one-touch of CAL key (internal weight model only).




Sophisticated Balance, the Professional's Choice

ViBRA LNA series always offers you the complete weighing solution. The capacity ranges from 620g to 31kg, the readability from 1mg to 0.1g. The advantages like quick response, clear

fluorescent display, tough housing, stylish design... ViBRA LNA series can be suitable for every occasions from laboratory, light&heavy industry, and jewelers.


Simple weighing



Carat for Jewelers



Percent




Parts Counting



Comparator



Density



Density measurement mode

To measure the density of the object is one of the most typical applications for the precision balance. ViBRA LNA series offers you the function to easily calculate the density from the measurement results. * the density measurement kit in the image is option.



Connection to the outside devices

ViBRA LNA series has RS232C as standard (two outputs) and can be easily connected to the printer, PC. You can keep the weighing results in the printed and/or electric forms.

