270021M21

HT(R) CEN / HT(R) CE Series
Specific Gravity Measurement Kit
Operation Manual

VIBRA

Preface

Thank you very much for having purchased Specific Measurement Kit for Electronic Analytical Balance HT(R) CEN / HT(R) CE series.

This document describes how to install and work with the Specific Measurement Kit.

Instructions

- The copyright of this document belongs to SHINKO DENSHI CO., LTD.. Reprinting or duplicating of all or part of this document without notice shall not be allowed.
- Please note that product improvement or modification may cause partial discrepancy between the product and the description of this document.
- The description of this document is subject to change without notice.
- This document has been created carefully. In case that, however, any error or imperfection is found by any chance, please let us know.
- Documents of which pages are missing or irregularly bound will be exchanged. Please inform the store where you purchased the product.
- Trouble related to the product or system will be dealt with in accordance with the individual maintenance contract. Please note, however, that we will not take responsibility for consequential trouble such as discontinuation of operation caused by the product trouble.
- **VIBRA** is the registered trademark of SHINKO DENSHI CO., LTD.. Company names and product names appearing in this document are the trademarks or registered trademarks of the respective company concerned.

Important Notice




- It should be known that this product contains potential danger. And so please be sure to observe this document when installing, operating or servicing this product.
- SHINKO DENSHI CO., LTD. will not take any responsibility for any injury or damage caused by the non-observance of this document or misuse or unauthorized modification of this product.

- Potential dangers are increasing in the industrial equipment industries due to the advent of new materials and processing methods, and speeding up of machines. It is impossible to foresee all situations related to these dangers. In addition, there are so many “impossible” and “don’ts” and so writing all of them in the operation manual is impossible. Therefore, it is safe to think that what is not written in the operation manual “cannot be performed” unless the operation manual positively writes “it is possible.” When performing installation, operation, maintenance or inspection of this product, not only observe what is written or indicated in this document or on the product surface but also pay adequate consideration to safety measures.
- The copyright of this document is held and reserved by SHINKO DENSHI CO., LTD.. Duplicating or disclosing its drawings and engineering materials without prior approval of SHINKO DENSHI CO., LTD. in writing is not permitted.
- For any question or further information concerning this document, please contact the store where you purchased the product with its model (type) name informed.
- Manufacturer: SHINKO DENSHI CO., LTD.
Address: 1-52-1 Itabashi, Itabashi-ku, Tokyo 173-0004 Japan

How to use this document

■ Symbols used in this document

Understand the meaning of the following symbols and observe the instructions of this document.

Symbols	Meaning
	Used for high risk point concerning the operations that may lead to death or severe physical injury to persons if proper precautions are not taken.
	Used for warning concerning the operations that may lead to death or severe physical injury to persons, if proper precautions are not taken.
	Used for caution concerning operations that may lead to a light physical injury to persons or damage of the products/facilities if proper precautions are not taken.
	Used for notation for avoiding from deletion, overwrite of the weighing data or for accurate weighing and appropriate usage of the equipment.
	Used for referenced information which is useful for product operation.
	Used for "Prohibition" items
	Used for "Mandatory" items requiring positive action
	This symbol indicates the operation/specification in related to the verified balance for legal metrology.

This product/ The product/The balance	Refers to the product.
[On/Off] key	The name of an operation key located in front of the main unit is represented in square brackets "[]".
<message>	A message on the display is represented in angle brackets "< >".
<<F1>>	"Free key" or "Shortcut" is represented in double angle brackets "<< >>".
Push the key	Signifies pushing lightly an operation key once.
Push the key long	Signifies keeping pushing an operation key until the designated indication appears.

Contents

Preface	i
Important Notice	ii
How to use this document	iii
Contents	iv
1 Prior to use	1
1-1 Operating precautions.....	1
1-2 For more accurate measurement.....	1
1-2-1 Precautions related to a specimen and medium liquid	1
1-2-5 Precautions related to the specific gravity measuring kit.....	1
1-3 Check for the articles contained in the box	2
2 Assembling and installation	3
3 Procedure for specific gravity measurement	5
3-1 Measuring specific gravity of the solid sample.....	6
3-1-1 Procedure for measuring the solid specimen	6
3-1-2 Switching the display at solid specific gravity mode	9
3-2 Measuring specific gravity of the liquid sample	10
3-2-1 Procedure for measuring specific gravity of the liquid sample.....	10
3-2-2 Switching the display at liquid specific gravity mode	14
4 Specifications	15
5 Print Details	16

1 Prior to use

In addition to the operating precaution of the balance body, following precautions shall be taken care for the specific gravity measurement kit.

1-1 Operating precautions


■ **Do not use/measure flammable liquids which generate explosive gases, vapor or fog.**

That may cause explosion or fire. A balance with explosion-proof structure shall be used to handle with such a flammable liquid.

■ **Do not install the balance on unstable or vibrating location, do not move the balance with a water tank set on it, especially when handling with hazardous liquids.**

Specific gravity measurement kit is not fixed to the balance, so that improper installation, carrying or operation may cause spilling of the hazardous liquids and lead to serious accidents.

■ **Take care not to spill water or liquid in the windshield.**

If a large amount of water or liquid is spilled in the windshield, they may enter inside the balance and cause shortsircuit or failure.

1-2 For more accurate measurement

1-2-1 Precautions related to a specimen and medium liquid

Unsuitable specimen	→ Specimen that is pneumatic, foamable or solubility cannot be measured accurately with this product.
Magnetism	→ Specimens affected by strong magnetism repels diamagnetic liquid, and that lead to an error of the measuring result.
Foam/bubbles on the specimen	→ Foam or bubbles produced on the specimen or the glass weight causes an error of the measuring result. If water is used as the medium, add one or two drops of surface-active agent (liquid detergent for kitchen use). It will effectively prevent air bubbles from adhering to any measurement object.
Temperature difference	→ Difference in temperature between the specimen and the medium liquid generates convection flow within the water tank, causing a measurement error.
Accuracy of the water temperature	→ When water is used as the medium, the accuracy of temperature of the water input into the balance affects accuracy of the measurement. Check the water temperature as frequently as appropriate.





1-2-5 Precautions related to the specific gravity measuring kit

Undesired contact with measurement pan	→ You should be careful not to leave the thermometer in the water tank. It may hit the pan submerged in water to result in inaccurate measurements. When water temperature is to be measured, once remove the water tank from the balance.
--	--

1-3 Check for the articles contained in the box


The package box contains the following;

If anything missing or broken should be found, please inform the store where you purchased the product.

<p>① Water tank platform</p>  A white, octagonal platform with four small legs, used for supporting a water tank.	<p>② Measurement Pan</p>  A metal measurement pan assembly. It consists of a shallow upper pan on a stem, which is attached to a wire mesh underwater pan. Labels "Upper pan" and "Underwater pan" are present. <p>Upper pan</p> <p>Underwater pan</p>	<p>③ Pan cradle</p>  A metal pan cradle with a vertical support and a horizontal arm that holds a shallow pan.
<p>④ Height adjuster</p>  A white L-shaped metal height adjuster and a small white plastic cap with a brass pin.	<p>⑤ Water tank</p>  A clear glass water tank.	<p>⑥ Glass weight</p>  A small glass weight with a circular base and a central knob.
<p>⑦ Tweezers</p>  A pair of metal tweezers.	<p>⑧ Thermometer</p>  A yellow thermometer and a silver thermometer.	<p>⑨ Operation manual</p>  A white rectangular operation manual.

2 Assembling and installation

- 1** Remove the weighing pan and pan base of the balance.

- 2** Attach the height adjuster to the pan cradle when needed.

- 3** Attach the pan cradle (turn the screw in the centre)

- 4** Attach the water tank platform.

5

Put the water tank in the position.

6

Set the measurement pan

3 Procedure for specific gravity measurement

Specific gravity is measured with “Specific gravity mode”. In the specific gravity mode, the ratio of the density of a substance to the density of water at its densest (4°C) for liquids is calculated.

Following modes for specific gravity measurement are available:

- Solid specific measurement mode: Measure solid sample with water or other medium liquid.
- Liquid specific measurement mode: Measure liquid sample with reference grass weight.

Legal
Metrology

Specific gravity mode is not legal for trade.

Note

- (1) The total weight of “Pan cradle”, “Measurement pan” and height adjuster is heavier than the weighing pan initially attached to the balance.
Therefore, the difference of the weight is deducted as tare value when the balance is used with this kit, and the maximum capacity to be measured in the balance decreases accordingly.
- (2) The diameter of the wire of the Measurement pan may affects the measurement result especially when volume or specific gravity of the solid sample is small.

Reference

Zero-point adjustment and tare-subtraction also can be executed by command input or external contact input. Please refer to “6 External input/output functions” of “HT/HTR CEN Series operation manual” or “HT/HTR CE Series operation manual”.

3-1 Measuring specific gravity of the solid sample


Outline of the procedure

1. Input the water temperature or the specific gravity of the medium liquid.
2. Compensate the buoyancy acting on the underwater part of the measurement pan by tare-subtraction.
3. Measure the specimen weight in the air.
4. Measure the specimen weight in the water/liquid.
5. The specific gravity of the specimen is displayed.

3-1-1 Procedure for measuring the solid specimen

1

Select the specific gravity mode.

Push [Menu] key, then push [F1]-[F4] keys to go to <11 MODE>

Push [F4] key to change the setting value.

Push [F1]/[F2] key to select.


SPGR: specific gravity mode

Push [F4] key to fix.

Push [Menu] key to shift to the specific gravity mode.

2

Select the measurement mode.

Push [F3]/[F4] key to select the measurement mode.

SOLID: Solid specific gravity mode

LIQ: Liquid specific gravity mode

3

Select the medium liquid.


Push [F3]/[F4] key to select the medium liquid.

OTHER: Liquid other than water having a known specific gravity value

H2O: water

4

Run medium liquid to the water tank and take the temperature.

Run medium liquid such as water to the attached water tank.


When water is selected as medium liquid, take the temperature and record it.

When other liquid is selected, take the temperature and look up the specific gravity at that temperature.

5


Input the specific gravity of the medium liquid or the temperature of the water.

<OTHER>: Liquid other than water

Enter the specific gravity of the medium liquid and push [F4] key to fix.

<H2O>: Water

Enter the temperature of the water and push [F4] key to fix.

When the specimen sinks


6 Assemble the specific gravity measurement kit.

Refer to "2 Assembling and installation".

7 Measure the specimen weight in the air.


Press [Zero/Tare] key to execute tare subtraction.

Load the specimen on the upper pan to measure the weight of the sample in the air, then push [F4] key to record it.

8 Measure the specimen weight in the water/liquid.

Put the specimen on the underwater pan, then push [F4] key to record it.

9

The specific gravity of the specimen (for the 4°C water) is automatically calculated and displayed.

The result can be output by pressing [Output] key.


When the specimen floats

6 Assemble the specific gravity measurement kit.

Refer to "2 Assembling and installation". When the specific gravity of the specimen is smaller than that of the medium liquid (the specimen floats), use the height adjuster.


7 Place a weight on the upper pan and execute tare subtraction.

Place a weight on the upper pan. Push [Zero/Tare] key to execute tare subtraction.


8 Measure the specimen weight in the air.

Load the specimen on the upper pan to measure the weight of the sample in the air, then push [F4] key to record it.

9 Measure the specimen weight in the water/liquid.

Put the specimen UNDER the underwater pan, leaving the weight on the upper pan. Then push [F4] key to record it.

10


The specific gravity of the specimen (for the 4°C water) is automatically calculated and displayed. The result can be output by pressing [Output] key.

Note

Be sure that the measurement pan does not come off from the pan cradle.

3-1-2 Switching the display at solid specific gravity mode

1 Push [F1]-[F4] keys to switch the display.


3-2 Measuring specific gravity of the liquid sample

Specific gravity of liquid is measured by using the accessory glass weight as a reference sample. In this measurement procedure, specific gravity of the glass weight is measured before measuring the specific gravity of liquid. The specific gravity value of the glass weight will be saved once it is measured. So, the corresponding operation may be skipped in subsequent measurements. You can also use other reference samples in place of the accessory glass weight. In doing so, you must perform a series of measurement steps to first measure the specific gravity of the reference sample and then to measure liquid specific gravity. Subsequently, you can skip the specific gravity measurement for the sample. Note that using another sample will overwrite the specific gravity value of the glass weight. If you are to use the glass weight again, you must measure the specific gravity of the glass weight before measuring liquid specific gravity.

3-2-1 Procedure for measuring specific gravity of the liquid sample

1

Select the specific gravity mode.

Push [Menu] key, then push [F1]-[F4] keys to go to <11 MODE>

Push [F4] key to change the setting value.

Push [F1]/[F2] key to select.


SPGR: specific gravity mode

Push [F4] key to fix.

Push [Menu] key to shift to the specific gravity mode.

2

Select the measurement mode.


Push [F3]/[F4] key to select the measurement mode.

SOLID: Solid specific gravity mode

LIQ: Liquid specific gravity mode

3

Select execution of sampling of reference weight.


Press the [F3]/[F4] key.

SKIP: Use stored data or input numerical value by key operation.

NO: Execute the sampling of the reference weight

When <SKIP> is selected:

4 Input the weight of the reference sample


➔ Numeric value input ➔

The weight of the last used reference sample is displayed.
When new reference sample to be used, input the weight of the reference sample.
Press the [F4] key.

5 Input the specific gravity of the reference sample

➔ Numeric value input ➔

The specific gravity of the last used reference sample is displayed.
When new reference sample to be used, input the specific gravity of the sample.
Press the [F4] key.
<CHANGE LIQUID> is displayed for several seconds.

6 Run the liquid specimen to the water tank.


7 Assemble the specific gravity measurement kit.

Refer to "2 Assembling and installation".


8 Submerge the underwater pan and execute tare subtraction.

The underwater pan into the liquid specimen and push [Zero/Tare] key to execute tare subtraction.

9 Measure the weight of the reference sample in the liquid specimen.

Put the reference sample on the pan submerged in the liquid specimen.
Then push [F4] key to record it.


10

The specific gravity of the liquid specimen (for the 4°C water) is automatically calculated and displayed.
The result can be output by pressing [Output] key.

When <NO> is selected:


4 Select medium liquid for measuring the reference sample.

Push [F3]/[F4] key to select the medium liquid.

OTHER: Liquid other than water having a known specific gravity value
H2O: water

5 Run medium liquid to the water tank and take the temperature.


Run medium liquid such as water to the attached water tank.

When water is selected as medium liquid, take the temperature and record it.

When other liquid is selected, take the temperature and look up the specific gravity at that temperature.


6 Input the specific gravity of the medium liquid or the temperature of the water.

<OTHER>: Liquid other than water

Enter the specific gravity of the medium liquid and push [F4] key to fix.

<H2O>: Water

Enter the temperature of the water and push [F4] key to fix.


7 Assemble the specific gravity measurement kit.

Refer to "2 Assembling and installation".


8 Measure weight of the reference sample in the air.

Press [Zero/Tare] key to execute tare subtraction.

Load the reference sample on the upper pan to measure the weight of the sample in the air, then push [F4] key to record it.

9 Measure weight of the reference sample in the water/liquid.

Put the reference sample on the underwater pan, then push [F4] key to record it.
 <CHANGE LIQUID> is displayed for several seconds.

10 Replace the liquid in the tank.

Remove the reference sample, measurement pan and liquid in the tank, then dry them.


11 Run the liquid specimen to the water tank.

12 Assemble the specific gravity measurement kit.

Refer to "2 Assembling and installation".


13 Submerge the underwater pan and execute tare subtraction.

The underwater pan into the liquid specimen and push [Zero/Tare] key to execute tare subtraction.

14 Measure the weight of the reference sample in the liquid specimen.


Put the reference sample on the pan submerged in the liquid specimen.
 Then push [F4] key to record it.

The specific gravity of the liquid specimen (for the 4°C water) is automatically calculated and displayed. The result can be output by pressing [Output] key.

3-2-2 Switching the display at liquid specific gravity mode

1 Push [F1]-[F4] keys to switch the display.

4 Specifications

<p>·Solid specific gravity measurement mode</p>
<p>1. Media liquid settings</p> <ul style="list-style-type: none">- Water: Setting Value: Water temperature Displayed Unit: °C Setting Range: 0.0 to 99.9 °C- Other than water: Setting value: Specific gravity of the media liquid Setting Range: Minimum 0.0001 * Result of liquid specific gravity measurement can be used. (Refer to the operation procedure of "3-2 Measuring specific gravity of the liquid sample")
<p>2. Readability of specific gravity result</p> <p>d = 0.0001</p>
<p>·Liquid specific gravity measurement mode</p>
<p>1. Media liquid settings for determining the specific gravity of the reference sample:</p> <ul style="list-style-type: none">- Water: Setting value: Temperature of the water used Same as with solid specific gravity- Other than water: Setting value: Specific gravity of the liquid used Same as with solid specific gravity
<p>2. Reference sample settings:</p> <p>Setting Value: Mass of glass weight or other reference sample [g] Setting Value: Specific gravity of glass weight or other reference sample</p> <ul style="list-style-type: none">- Setting Range: Minimum 0.0001 <p>* Once a series of measurements, i.e., measurement of the specific gravity of glass weight followed by the measurement of liquid specific gravity, is performed, the necessary settings are automatically recorded, and they are reusable in later works. It is also allowed to set the parameters manually. In this case, the steps to measure the specific gravity of the glass weight is skipped. When the specific gravity of the glass weight is measured, the parameters already set are overwritten by the new measurement. If any other sample is used instead of the glass weight, the mass and specific gravity values of the alternative sample are set as the corresponding parameter values.</p>
<p>3. Readability of specific gravity result:</p> <p>d = 0.0001</p>

5 Print Details

[Print samples for CSP-240]

Mode: Solid

Medium: Water

Output: Specific gravity of
the sample / Weight
of the sample /
Temperature of the
water

SAMPLE SP GR
SAMPLE WEIGHT
WATER TEMP

Mode: Solid

Medium: Other than water

Output: Secific gravity of the
sample / Weight of
the sample / Specific
gravity of the media
liquid

SAMPLE SP GR
SAMPLE WEIGHT
MED. LIQ SP GR

Mode: Liquid

Medium: NA

Output: Liquid specific
gravity of the
sample

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