


AB Series

Tuning-fork Precision balance

High Precision in
Affordable Budget



JQA-2834
TSUKUBA



ISO 14001:2004
ER-00130

Specifications

Model	AB 323	AB 623	AB 1202	AB 3202	AB 12001
Capacity	320g	620g	1200g	3200g	12000g
Readability (d)	0.001g	0.001g	0.01g	0.01g	0.1g
Repeatability (sd.)	0.001g	0.001g	0.01g	0.01g	0.1g
Non-Linearity (typ.)	±0.003g	±0.004g	±0.02g	±0.03g	±0.2g
Stabilization time (typ.)	3.0 sec.	3.5 sec.	3.0 sec.	3.5 sec.	3.0 sec.
Pan size	140mm dia		190×190mm		
Calibration			External		
Output			RS232C		
Dimension	293×202×206mm *1		293×196×89mm		
Weight	approx.3.5kg		approx.2.6kg		

*1 including windshield

